



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Julabo 1.5 m metal hose, triple insulated, 2 connections M24x1.5 inside (-100...+350

€612.00
plus VAT &
Shipping

Product Images



Description

Technical data: - Working temperature range: -100°C ... +350°C - Material: stainless steel 1.4404 - Insulation: special triple insulation. The ends are covered with silicone flaps. - Connection: M24x1.5 thread - Length: 1.5 m - Outer diameter: 45 mm - Clearance: 12mm - Maximum pressure at 20°C: 6 bar - Maximum pressure at 350°C: 4 bar - Min. bending radius: 220mm - Fire rating: UL94-HB, UL94-V2

Additional Information

No.	JN-0156
Manufacturer (Brand)	Julabo
EAN	4058072844561
Transport temperature	Room temperature
Net weight	2.20000 kg
Country of origin	Switzerland
Customs tariff number	83071000

