



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lauda metal hose M30X 100S M30x1,5-M30x1,5
DN 20 NL=1000mm / Temp.-ber. -100...350**

**€706.00
plus VAT &
Shipping**

Product Images

Description

Metal hose M30X 100S M30x1,5-M30x1,5 DN 20 NL=1000mm / Temp.-ber. -100...350°C Performance features: Corrugated metal hose for temperature control with different temperature control fluids Multi-layer insulation, suitable for large temperature ranges Identification of the hose with permissible temperature range and maximum operating pressure Technical data: Operating temperature range: -100 ... 350 °C Operating pressure: -0.3 ... 10.0 bar in operating temperature range Hydraulic connection 1: M30 x 1.5 Hydraulic connection 2: M30 x 1.5 Length: 100 cm Inner diameter: 20 mm Outer diameter: 76 mm Material inside: Stainless steel Material outside: Smooth silicone hose with silicone end caps Vacuum resistance: 100 % for 0 to 50 °C Classification: UL 94 HB

Additional Information

No.	LA-0852
Manufacturer (Brand)	Lauda
EAN	4058072882686
Gross weight	5.5
Transport temperature	Room temperature
Net weight	3.50000 kg
Country of origin	Switzerland
Packing width	1.00000 m
Packing height	0.07500 m
Packing depth	0.07500 m
Customs tariff number	83071000

