



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lauda metal hose M38X 200S M38x1,5-M38x1,5
DN 25 NL=2000mm / Temp.-ber. -100...350**

**€1,270.00
plus VAT &
Shipping**

Product Images

Description

Metal hose M38X 200S M38x1.5-M38x1.5 DN 25 NL=2000mm / Temp.-ber. -100...350°C Performance features: Corrugated metal hose for temperature control with different temperature control fluids Multi-layer insulation, suitable for large temperature ranges Identification of the hose with permissible temperature range and maximum operating pressure Technical data: Operating temperature range: -100 ... 350 °C Operating pressure: -0.3 ... 10.0 bar in operating temperature range Hydraulic connection 1: M38 x 1.5 Hydraulic connection 2: M38 x 1.5 Length: 200 cm Inner diameter: 25 mm Outer diameter: 78 mm Material inside: Stainless steel Material outside: Smooth silicone hose with silicone end caps Vacuum resistance: 100 % for 0 to 50 °C Classification: UL 94 HB

Additional Information

No.	LA-0856
Manufacturer (Brand)	Lauda
EAN	4058072882723
Transport temperature	Room temperature
Net weight	6.32000 kg
Country of origin	Switzerland
Packing width	2.00000 m
Packing height	0.08000 m
Packing depth	0.08000 m
Customs tariff number	83071000

