



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Lauda metal hose M38X 300S M38x1,5-M38x1,5  
DN 25 NL=3000mm / Temp.-ber. -100...350**

**€1,760.00  
plus VAT &  
Shipping**

## Product Images

---

## Description

---

Metal hose M38X 300S M38x1.5-M38x1.5 DN 25 NL=3000mm / Temp.-ber. -100...350°C Performance features: Corrugated metal hose for temperature control with different temperature control fluids Multi-layer insulation, suitable for large temperature ranges Identification of the hose with permissible temperature range and maximum operating pressure Technical data: Operating temperature range: -100 ... 350 °C Operating pressure: -0.3 ... 10.0 bar in operating temperature range Hydraulic connection 1: M38 x 1.5 Hydraulic connection 2: M38 x 1.5 Length: 300 cm Inner diameter: 25 mm Outer diameter: 78 mm Material inside: Stainless steel Material outside: Smooth silicone hose with silicone end caps Vacuum resistance: 100 % for 0 to 50 °C Classification: UL 94 HB

## Additional Information

---

No.	LA-0857
Manufacturer (Brand)	Lauda
EAN	4058072882730
Transport temperature	Room temperature
Net weight	8.50000 kg
Country of origin	Switzerland
Customs tariff number	83071000

