



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Double connector T-shape f. tube I.-Ø 1.5 mm, 10 pcs./pack

**€25.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Double-T connector made of polypropylene, autoclavable up to 121°C.

The cone has a special shape so that hoses cannot slip down.

Also suitable for high-pressure applications.

In a bag of 10 pcs

## Additional Information

---

No.	2-0997
Manufacturer (Brand)	neoLab
EAN	4058072048891
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Form	Double T
Connector properties	Miniature version
for DM hose inside 1	1.5 mm
Symmetry	for same hose DM
Type hose connector	Double T-connector

