



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Double angle connector f. tube inner Ø 2.25 mm, 10 pcs./pack

Special Price

€15.00 was

€22.50

**plus VAT &
Shipping**

Product Images



Description

Double angle connector made of polypropylene, autoclavable up to 121°C.

The cone has a special shape so that hoses cannot slip down.

Also suitable for high-pressure applications.

In a bag of 10 pcs

Additional Information

No.	2-0990
Manufacturer (Brand)	neoLab
EAN	4058072048822
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Form	Double angle
Symmetry	for same hose DM
Type hose connector	Double angle connector
for DM hose inside 1	2.25 mm

