



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® 3-fold connector f. tube inner Ø 2.25 mm, 10 pcs./pack


Special Price

**€15.00** was

**€22.50**

**plus VAT &  
Shipping**

### Product Images



## Description

Hose connector, 3-fold, made of polypropylene, autoclavable up to 121°C.

The cone has a special shape so that hoses cannot slip down.

Also suitable for high-pressure applications.

In a bag of 10 pcs

## Additional Information

No.	2-0984
Manufacturer (Brand)	neoLab
EAN	4058072048761
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Form	straight angled
for DM hose inside 1	2.25 mm
Symmetry	for same hose DM
Type hose connector	Hose connector (3-fold)

