



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Connector Luer-Cone male, olive 1.6-2.6 mm, PP, 10 pcs/pack

€26.10
plus VAT &
Shipping

Product Images



Description

Tube Luer connector with olive/luer cone (m)

Made of natural colored PP, slightly transparent, **autoclavable**. Outstanding quality with constant inner diameter, causes uniform flow rates and tight connection. The olives are worked on the outside of the sealing cone in such a way that solid tubing can also be attached. They expand the hose up to 1.5 times its inner \varnothing .

In a bag of 10 pieces.

Additional Information

No.	2-1890
Manufacturer (Brand)	neoLab
EAN	4058072058678
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
for DM hose inside 1	1.6 to 2.6 mm
Gender Luer cone 1	male
Type Hose Luer Connector	Olive/ Luer cone male

