



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Connector, LL female/ Olive 2.4-3.4 mm, PP, 10 pcs./pack

€26.10
plus VAT &
Shipping

Product Images



Description

Tube-Luer-Connector with olive and Luer-Lock-Cone (w)

Made of natural colored PP, slightly transparent, **autoclavable**. Excellent quality with constant inner diameter, causes uniform flow rates and tight connection. The olives are worked on the outside of the sealing cone in such a way that solid tubing can also be attached. They expand the hose up to 1.5 times its inner \varnothing .

In a bag of 10 pieces.

Additional Information

No.	2-1887
Manufacturer (Brand)	neoLab
EAN	4058072058647
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Gender Luer lock	female
Type Hose Luer Connector	Olive/ Luer lock female
for DM hose inside 1	2.4 to 3.4 mm

