



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® Connector, LL female/ Olive 4.0-5.0  
mm, PP, 10 pcs./pack**

**€26.10  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

### Tube-Luer-Connector with olive and Luer-Lock-Cone (w)

Made of natural colored PP, slightly transparent, **autoclavable**. Excellent quality with constant inner diameter, causes uniform flow rates and tight connection. The olives are worked on the outside of the sealing cone in such a way that solid tubing can also be attached. They expand the hose up to 1.5 times its inner  $\varnothing$ .

**In a bag of 10 pieces.**

## Additional Information

---

No.	2-1889
Manufacturer (Brand)	neoLab
EAN	4058072058661
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
for DM hose inside 1	4.0 to 5.0 mm
Gender Luer lock	female
Type Hose Luer Connector	Olive/ Luer lock female

