



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PA tube connector, f. Inner-Ø 3-5/10-13 mm, 10 pcs./pack

€6.03
plus VAT &
Shipping

Product Images



Description

Lightweight and unbreakable hose connectors

Non-toxic, suitable for use in medicine, biology and food testing.

Shape of seamless lumina ensures streamlined connection and flow. Temperature resistant up to 100°C

Packing unit: 10 pcs./bag

Additional Information

No.	2-0880
Manufacturer (Brand)	neoLab
EAN	4058072048334
autoclavable at 121°C	No
Color	transparent (milky)
Material	Polyamide (PA. non-specific)
Form	straight
for DM hose inside 1	3.0 to 5.0 mm
for DM hose inside 2	10.0 to 13.0 mm
Symmetry	for different hose DM
Type hose connector	Reducer (straight)

