



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® PA tube connector, f. Inner-Ø 3-5/7-10 mm, 10 pcs./pack

**€18.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Lightweight and unbreakable hose connectors

Non-toxic, suitable for use in medicine, biology and food testing.

Shape of seamless lumina ensures streamlined connection and flow. Temperature resistant up to 100°C

**Packing unit:** 10 pcs./bag

## Additional Information

---

No.	2-0881
Manufacturer (Brand)	neoLab
EAN	4058072048341
autoclavable at 121°C	No
Color	transparent (milky)
Material	Polyamide (PA. non-specific)
Form	straight
for DM hose inside 1	3.0 to 5.0 mm
for DM hose inside 2	7.0 to 10.0 mm
Symmetry	for different hose DM
Type hose connector	Reducer (straight)

