



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® PA tube connector f. inner Ø 7-10/10-13
mm, 10 pcs./pack**

**€18.20
plus VAT &
Shipping**

Product Images



Description

Lightweight and unbreakable hose connectors

Non-toxic, suitable for use in medicine, biology and food testing.

Shape of seamless lumina ensures streamlined connection and flow. Temperature resistant up to 100°C

Packing unit: 10 pcs./bag

Additional Information

No.	2-0884
Manufacturer (Brand)	neoLab
EAN	4058072048372
autoclavable at 121°C	No
Color	transparent (milky)
Material	Polyamide (PA. non-specific)
Form	straight
Symmetry	for different hose DM
Type hose connector	Reducer (straight)
for DM hose inside 1	7.0 to 10.0 mm
for DM hose inside 2	10.0 to 13.0 mm

