



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PA tube connector, Y-shape, f. Inner-Ø 4-6 mm, 10 pcs./pack

€9.10
plus VAT &
Shipping

Product Images



Description

Lightweight and unbreakable hose connectors

Non-toxic, suitable for use in medicine, biology and food testing.

Shape of seamless lumina ensures streamlined connection and flow. Temperature resistant up to 100°C

Packing unit: 10 pcs./bag

Additional Information

No.	2-0890
Manufacturer (Brand)	neoLab
EAN	4058072048396
autoclavable at 121°C	No
Color	transparent (milky)
Material	Polyamide (PA. non-specific)
Form	Y-shaped
for DM hose inside 1	4.0 to 6.0 mm
Symmetry	for same hose DM
Type hose connector	Y-connector

