



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Plastic tubing connector TPX, Y-shape, 2-5 mm, 10 pcs./pack

**€58.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

TPX resin is highly transparent, non-toxic and therefore can be used in food, medical and biological fields. Temperature resistance up to 200°C.

The constant inner diameter allows a uniform flow rate.

**PU:** 10 pcs./bag

## Additional Information

---

|                       |                              |
|-----------------------|------------------------------|
| No.                   | 2-1853                       |
| Manufacturer (Brand)  | neoLab                       |
| EAN                   | 4058072058302                |
| Transport temperature | Room temperature             |
| autoclavable at 121°C | Yes                          |
| Color                 | transparent (clear)          |
| Material              | Polymethylpentene (PMP. TPX) |
| Form                  | Y-shaped                     |
| for DM hose inside 1  | 2.0 to 5.0 mm                |
| Symmetry              | for same hose DM             |
| Type hose connector   | Y-connector                  |

