



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Plastic tubing connector TPX, Y-shape, 5-9 mm, 10 pcs./pack

€71.20
plus VAT &
Shipping

Product Images



Description

TPX resin is highly transparent, non-toxic and therefore can be used in food, medical and biological fields. Temperature resistance up to 200°C.

The constant inner diameter allows a uniform flow rate.

PU: 10 pcs./bag

Additional Information

No.	2-1854
Manufacturer (Brand)	neoLab
EAN	4058072058319
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylpentene (PMP. TPX)
Form	Y-shaped
for DM hose inside 1	5.0 to 9.0 mm
Symmetry	for same hose DM
Type hose connector	Y-connector

