



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Plastic tubing connector TPX, Y-shape, 8-15 mm, 10 pcs./pack

**€93.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

TPX resin is highly transparent, non-toxic and therefore can be used in food, medical and biological fields. Temperature resistance up to 200°C.

The constant inner diameter allows a uniform flow rate.

**PU:** 10 pcs./bag

## Additional Information

---

No.	2-1855
Manufacturer (Brand)	neoLab
EAN	4058072058326
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylpentene (PMP. TPX)
Form	Y-shaped
Symmetry	for same hose DM
Type hose connector	Y-connector
for DM hose inside 1	8.0 to 15.0 mm

