



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Plastic tubing connector TPX, transition piece, 5-8/8-15 mm, 10 pcs./pack

€66.50
plus VAT &
Shipping

Product Images



Description

TPX resin is highly transparent, non-toxic and therefore can be used in food, medical and biological fields. Temperature resistance up to 200°C.

The constant inner diameter allows a uniform flow rate.

PU: 10 pcs./bag

Additional Information

No.	2-1860
Manufacturer (Brand)	neoLab
EAN	4058072058371
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylpentene (PMP. TPX)
Form	straight
for DM hose inside 1	5.0 to 8.0 mm
for DM hose inside 2	8.0 to 15.0 mm
Symmetry	for different hose DM
Type hose connector	Reducer (straight)

