



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Connector straight, PP, olive on both sides for 2.4–3.4 mm, 10 pcs./pack

**€25.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Connector, straight, made of natural-colored PP**

**Both sides with the same olive size.** Slightly transparent, **autoclavable**. Excellent quality with constant inner diameter, resulting in uniform flow rates and tight connection. The olives are machined on the outside of the sealing cone so that rigid tubing can also be attached. They expand the hose up to 1.5 times its inner  $\varnothing$ .

**In a bag of 10 pieces.**

## Additional Information

---

No.	2-0959
Manufacturer (Brand)	neoLab
EAN	4058072048518
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Form	straight
for DM hose inside 1	2.4 to 3.4 mm
Symmetry	for same hose DM
Type hose connector	Hose connector (straight)

