



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Connector straight, PP, both sides for 4.8–5.8 mm, 10 pcs./pack

€25.50
plus VAT &
Shipping

Product Images



Description

Connector, straight, made of natural-colored PP

Both sides with the same olive size. Slightly transparent, **autoclavable**. Excellent quality with constant inner diameter, resulting in uniform flow rates and tight connection. The olives are machined on the outside of the sealing cone so that rigid tubing can also be attached. They expand the hose up to 1.5 times its inner \varnothing .

In a bag of 10 pieces.

Additional Information

No.	2-0967
Manufacturer (Brand)	neoLab
EAN	4058072048594
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Form	straight
for DM hose inside 1	4.8 to 5.8 mm
Symmetry	for same hose DM
Type hose connector	Hose connector (straight)

