



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® connector set, 72 parts in box

€192.00
plus VAT &
Shipping

Product Images



Description

A whole box full of good stuff! This 72-piece assortment has everything you need for everyday use.

The set includes the following parts:

- 15 straight connectors (PP)
- 6 tees (PP)
- 6 quick disconnect pieces (PE)
- 11 straight/stepped connectors (PP)
- 10 Y-pieces (PP)
- 24 polyethylene cup-type plugs (PE).

All parts are packed in a transparent polystyrene box with hinged lid.

Material PP: autoclavable and FDA grade

Additional Information

No.	6-2255
Manufacturer (Brand)	neoLab
EAN	4058072025595
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Number of parts in set	72
Type hose connector	Hose connector (straight) T-connector Reducer (straight) Y-connector

