



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Angle connector f. tube inner Ø 1.5
mm, 10 pcs./pack**

**€25.50
plus VAT &
Shipping**

Product Images



Description

Angle connector made of polypropylene, autoclavable up to 121°C.

The cone has a special shape so that hoses cannot slip down.

Also suitable for high-pressure applications.

In a bag of 10 pcs

Additional Information

| | |
|-----------------------|--------------------|
| No. | 2-0980 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072048723 |
| autoclavable at 121°C | Yes |
| Color | White |
| Material | Polypropylene (PP) |
| Form | angled |
| for DM hose inside 1 | 1.5 mm |
| Symmetry | for same hose DM |
| Type hose connector | Angle connector |

