



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Angle connector f. tube inner Ø 3 mm,  
10 pcs./pack**

**€25.50  
plus VAT &  
Shipping**

## Product Images



## Description

---

Angle connector made of polypropylene, autoclavable up to 121°C.

The cone has a special shape so that hoses cannot slip down.

Also suitable for high-pressure applications.

In a bag of 10 pcs

## Additional Information

---

No.	2-0982
Manufacturer (Brand)	neoLab
EAN	4058072048747
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Form	angled
Symmetry	for same hose DM
Type hose connector	Angle connector
for DM hose inside 1	3.0 mm

