



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Y-tube connector f. 1.6-2.6 mm, PP,  
autoclavable, 10 pcs/pack**

**€25.50  
plus VAT &  
Shipping**

## Product Images



## Description

---

Y-connector made of natural-colored PP

Same olive size on all sides. Slightly transparent, autoclavable. Outstanding quality with constant inner diameter, causes uniform flow rates and tight connection. The olives are machined on the outside of the sealing cone so that solid tubing can be attached. They expand the hose up to 1.5 times its inner  $\varnothing$ .

In a bag of 10 pieces.

## Additional Information

---

No.	2-0956
Manufacturer (Brand)	neoLab
EAN	4058072048488
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Form	Y-shaped
for DM hose inside 1	1.6 to 2.6 mm
Symmetry	for same hose DM
Type hose connector	Y-connector

