



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Y-tube connector f. 3.2-4.2 mm, PP,
autoclavable, 10 pcs/pack**

**€25.50
plus VAT &
Shipping**

Product Images



Description

Y-connector made of natural-colored PP

Same olive size on all sides. Slightly transparent, autoclavable. Outstanding quality with constant inner diameter, causes uniform flow rates and tight connection. The olives are machined on the outside of the sealing cone so that solid tubing can be attached. They expand the hose up to 1.5 times its inner \varnothing .

In a bag of 10 pieces.

Additional Information

No.	2-0958
Manufacturer (Brand)	neoLab
EAN	4058072048501
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Form	Y-shaped
Symmetry	for same hose DM
Type hose connector	Y-connector
for DM hose inside 1	3.2 to 4.2 mm

