



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® "Wendoflex" electrode holder

Special Price

€44.90 was

€62.90

**plus VAT &
Shipping**

Product Images



Description

Electrode holder "Wendoflex"

Inexpensive electrode holder with flexible metal hose and plastic clamps for electrodes up to $d = 10$ mm. 3 holders. The flexible hose allows the electrode to be placed in any desired position.

Total length: 380 mm

Additional Information

No.	2-4150
Manufacturer (Brand)	neolabLine
EAN	4058072002275
Transport temperature	Room temperature
suitable for	Electrodes Thermometer
Number of holders	1
Length rod / arm	380 mm
Material rod / arm	Steel chrome plated
Tripod type	Tripod with flexible arm
with holders for	Electrodes Thermometer

