



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Tube holder for suction devices, with magnet

**€61.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

This **magnetic hose holder** was developed in cooperation with the University Hospital of Tübingen and can be attached to the side wall of the sterile workbench, for example.

Because hygiene is the be-all and end-all when working under the workbench. With the counter magnet included in the scope of delivery, the hose holder can even be attached to glass panes.

Ideal e.g. for the suction hose of the suction systems order no. 2-9335/36 or 2-9339.

## Additional Information

---

No.	7-3047
Manufacturer (Brand)	neoLabLine
EAN	4058072038564
Transport temperature	Room temperature
Color	Black
Material	Steel powder coated
suitable for	Hoses
Type mounting	magnetic

