



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz holder for grid stands, complete with
tension lock and supports, grid stand grid 0.707**

**€1,436.14
plus VAT &
Shipping**

Product Images



Description

Holder for grid stands

For mounting attachments on reaction vessels and for hanging in grid stands, complete with clamping lock and holders, material: stainless steel (clamping lock), steel, powder-coated (holder)

Additional Information

No.	LZ-1887
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3890
EAN	4058072316259
Transport temperature	Room temperature
Color	dark gray
Material	Steel powder coated stainless steel
suitable for	Stands for reaction apparatus (LENZ) Vessels with LF 150
MAX span	707 mm
Span MIN	707 mm
Type sleeve	Bracket for grid tripods
Type connection	parallel

