



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® stainless steel quick release fasteners, with 3 retaining segments, suitable for flange DN 100

€83.13
plus VAT &
Shipping

Product Images

Description

Quick-release closure, DN 100 The quick-release closure closes flat-flange vessels and secures the lid. Slipping or uncontrolled opening is not possible. The closure is made of high-quality stainless steel and has three retaining segments. Special features and advantages: With retaining segments Fast closing due to quick-release closure Application: Quick-release ring for closure, accessories for plan flange vessels

Additional Information

No.	SD-2590
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	132907146
EAN	4058072362263
Manufacturer EAN	4032051061446
Gross weight	0.24
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	0.21000 kg
Material	Stainless steel
sterile	No
Country of origin	Great Britain
Customs tariff number	73269098

