



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® stainless steel quick release fasteners, with 3 retaining segments, suitable for flange DN 120

€87.79
**plus VAT &
Shipping**

Product Images



Description

Quick-release closure, DN 120 The quick-release closure closes flat-flange vessels and secures the lid. Slipping or uncontrolled opening is not possible. The closure is made of high-quality stainless steel and features three retaining segments. Special features and advantages: With retaining segments Fast closing due to quick-release closure
Application: Quick-release ring for closure, accessories for plan flange vessels

Additional Information

No.	SD-2591
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	132907151
EAN	4058072362270
Manufacturer EAN	4032051061460
Gross weight	0.244
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	0.24400 kg
Material	Stainless steel
sterile	No
Country of origin	Great Britain
Packaging volume	0.00059200 ccm
Customs tariff number	73269085

