




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Lifting Platform Poly-Jaque

€431.50
plus VAT &
Shipping

Product Images



Description

Lightweight but stable lifting platform, corrosion and heat resistant, chemical resistant.

Made of glass fiber reinforced polymer, loadable up to 7 kg.

Continuously adjustable in height from 90 to 305 mm via an adjustment screw.

PTFE-coated pivot screw.

Storage surface: 165 x 159 mm

Additional Information

No.	6-2565
Manufacturer (Brand)	neoLab
EAN	4058072026301
Transport temperature	Room temperature
Material	Polycarbonate (PC)
Wide	165 mm
Length	159 mm
Lifting platform features	with special chemical resistance
Height [MAX]	305 mm
Height [MIN]	90 mm
loadable up to	7 kg

