



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® spring clips 25 mm wide, 6 pcs./pack

€7.63
plus VAT &
Shipping

Product Images



Description

For many things that need to be held quickly and easily.

The spring clips are very strong, but can be easily bent open with the bracket.

The eyelet-shaped shackle also makes these clips suitable for hanging papers.

25 mm wide: locking width 9 mm

Additional Information

No.	2-4027
Manufacturer (Brand)	neoLab
EAN	4058072001797
Transport temperature	Room temperature
Color	black silver
Material	Steel powder coated
Wide	25 mm
suitable for	variable applications
Material clamp	Steel powder coated
MAX span	9 mm
Span MIN	0 mm

