



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® spring steel clamp, VA 18/8, for 1.5 ml, 10 mm

**€6.19** As low as  
**€5.60**  
plus VAT &  
Shipping

### Product Images



## Description

---

Screw-on clamps for a variety of purposes, for fastening small parts, hoses, etc.. The clamp that does not break or crack, does not rust and is ball polished. With hole for screwing.

## Additional Information

---

No.	1-2131
Manufacturer (Brand)	neoLab
EAN	4058072034214
Color	silver
Material	Stainless steel
suitable for	Centrifuge tubes tubes
Material clamp	Stainless steel
MAX span	10 mm
for clamping to	Turntable for rotators

