



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® spring steel clamp, galvanized steel, for 50 ml, 30 mm

€8.33 As low as

€7.27

**plus VAT &
Shipping**

Product Images



Description

Screw-on clamps for a variety of purposes, for fastening small parts, hoses, etc.. The clamp that does not break or crack, does not rust and is ball polished. With hole for screwing.

Additional Information

No.	1-2133
Manufacturer (Brand)	neoLab
EAN	4058072034238
Color	silver
Material	Steel galvanized
suitable for	Centrifuge tubes tubes
MAX span	30 mm
Material clamp	Steel galvanized
for clamping to	Turntable for rotators

