



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Corning clamps for 250 mL Erlenmeyer flasks

€45.92
plus VAT &
Shipping

Product Images



Description

Stainless steel clamps suitable for Axygen® incubator shakers (49 and 71 liters); volume: 250 ml

Additional Information

No.	CQ-0077
Manufacturer (Brand)	Corning
old neoLab article no.	C-2083
EAN	4058072225476
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
suitable for	Erlenmeyer flask
Material clamp	Stainless steel
for vessels with volume up to	250 ml
for clamping to	Shaker platform

