



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Corning clamps for 1 L Erlenmeyer flask

€45.92
plus VAT &
Shipping

Product Images



Description

Stainless steel clamps suitable for Axygen® incubator shakers (49 and 71 liters); volume: 1000 ml

Additional Information

| | |
|-------------------------------|------------------|
| No. | CQ-0079 |
| Manufacturer (Brand) | Corning |
| old neoLab article no. | C-2085 |
| EAN | 4058072225490 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Stainless steel |
| suitable for | Erlenmeyer flask |
| for vessels with volume up to | 1000 ml |
| Material clamp | Stainless steel |
| for clamping to | Shaker platform |

