



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Corning rack holder for 50ml test tubes

**€50.40**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Holder for 50 ml tubes

## Additional Information

---

No.	CQ-0134
Manufacturer (Brand)	Corning
old neoLab article no.	C-2079
EAN	4058072225735
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
suitable for	Centrifuge tube
for vessels with volume up to	50 ml
Material clamp	Stainless steel
for model shaker/ rotator	Axygen incubator shaker 71 liters (C-2068)
for clamping to	Shaker platform

