



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Corning plastic clamp, 500ml

€51.63
plus VAT &
Shipping

Product Images



Description

Plastic clamps suitable for Axygen® incubator shaker Benchtop; 500 ml

Additional Information

No.	CQ-0116
Manufacturer (Brand)	Corning
old neoLab article no.	C-2059
EAN	4058072225568
Transport temperature	Room temperature
Color	Black
Material	Acrylonitrile butadiene styrene (ABS)
suitable for	Erlenmeyer flask
for vessels with volume up to	500 ml
Material clamp	Acrylonitrile butadiene styrene (ABS)
for clamping to	Shaker platform

