



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® Clamp for 50 ml

€12.00
plus VAT &
Shipping

Product Images



Description

Clamps for fixing Erlenmeyer flasks of various volumes on the universal platform of the Sunlab orbital shaker.

Order No. D-8023: for 50 ml

Additional Information

No.	D-8023
Manufacturer (Brand)	sunlab
EAN	4058072071011
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
suitable for	Erlenmeyer flask
for vessels with volume up to	50 ml
Material clamp	Stainless steel
for model shaker/ rotator	Sunlab Orbital Shaker (D-8020)
for clamping to	Shaker platform

