



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## sunlab® Clamps for 1.5 ml centrifuge tubes suitable for Rotator D-8010/ D-8101

Special Price

**€8.00** was

**€39.90**

**plus VAT &  
Shipping**

### Product Images



## Description

---

### Sunlab clamps

for attaching various vessels to the rotator.

### Order No. D-8011

For 1.5 ml reaction tubes, suitable for rotator D-8010/D-8100

## Additional Information

---

No.	D-8011
Manufacturer (Brand)	sunlab
EAN	4058072017965
Transport temperature	Room temperature
Color	White
Material	Polyamide (PA. non-specific)
suitable for	Reaction vessels
MAX span	8 mm
Span MIN	12 mm
for vessels with volume up to	1.5 ml
Material clamp	Polyamide (PA. non-specific)
for model shaker/ rotator	Sunlab Rotator with turntable, analog (D-8010) Sunlab Rotator variable equipment (D-8100)
for clamping to	Turntable for rotators

