



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® Clamps for 15 ml centrifuge tubes suitable for Rotator D-8010/ D-8101

Special Price

€8.00 was

€39.90

**plus VAT &
Shipping**

Product Images



Description

Sunlab clamps

for attaching various vessels to the rotator.

Order No. D-8012

For 15 ml centrifuge tubes, suitable for rotator D-8010/D-8100

Additional Information

No.	D-8012
Manufacturer (Brand)	sunlab
EAN	4058072017972
Transport temperature	Room temperature
Color	White
Material	Polyamide (PA. non-specific)
suitable for	Reaction vessels Centrifuge tubes
MAX span	15 mm
Span MIN	20 mm
for vessels with volume up to	15 ml
Material clamp	Polyamide (PA. non-specific)
for model shaker/ rotator	Sunlab Rotator with turntable. analog (D-8010) Sunlab Rotator variable equipment (D-8100)
for clamping to	Turntable for rotators

