




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz ground joint clamp (POM) for NS 45

€38.40
plus VAT &
Shipping

Product Images



Description

Ground joint clamp made of POM.
Individual contact pressure due to stable locking linkage.
Temperature resistant up to +150°C.

Additional Information

No.	LZ-0138
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2137
EAN	4058072298807
Transport temperature	Room temperature
Material	Polyoxymethylene (POM)
TBST MAX	150 °C
for NS	45/40

