



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell quick holder for vessel walls

€48.05 As low
as
€36.69
plus VAT &
Shipping

Product Images



Description

For thermometers, test tubes with samples and the like,
for easy clamping to beakers, Erlenmeyer flasks, glass tubs, etc.
With five different insertion options.
All metal parts made of V2A, sturdy workmanship.

Additional Information

No.	LL-1491
Manufacturer (Brand)	Amarell
old neoLab article no.	1-7065
EAN	4058072288273
Color	silver
Material	Stainless steel
suitable for	Test tubes Thermometer
MAX span	3 mm
Span MIN	0 mm
Material clamp	Stainless steel
for clamping to	Beakers Pistons

