



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz clamping lock for mounting attachments on reaction vessels, for laboratory flange LF 60

€73.72
plus VAT &
Shipping

Product Images



Description

Clamping caps

For mounting attachments on reaction vessels, for flat flanges LF 60

Additional Information

No.	LZ-1873
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3876
EAN	4058072316112
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Form	around
suitable for	Vessels with LF 60
MAX span	60 mm
Span MIN	60 mm

