



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz clamping lock with stand rod, for laboratory flange LF 150, lock Polyester, glass fiber reinforced

€978.24
plus VAT &
Shipping

Product Images



Description

Clamps with stand rod

For fastening attachments to reaction vessels and mounting on a stand, material: polyester, glass fiber reinforced (closure), steel (rod)

Additional Information

| | |
|------------------------|---------------------|
| No. | LZ-1879 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-3882 |
| EAN | 4058072316174 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Stainless steel |
| Form | around |
| Length | 300 mm |
| suitable for | Vessels with LF 150 |
| MAX span | 150 mm |
| Span MIN | 150 mm |
| for DM Tripod Rod MAX | 16 mm |
| for DM tripod rod MIN | 10 mm |
| for clamping to | Tripod rods |

