



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Thermometer clips, 3 pcs./pack

€23.20
plus VAT &
Shipping

Product Images



Description

Simply place the thermometer clip on the rim of a beaker or other vessel. It clamps onto the rim and holds the thermometer vertically, close to the vessel wall. Slotted clamping device for thermometers with different diameters.

Material: Polycarbonate

Temperature resistance: up to 135°C.

Additional Information

No.	1-3090
Manufacturer (Brand)	neoLab
EAN	4058072035006
Color	transparent (clear)
Material	Polypropylene (PP)
suitable for	Thermometer
MAX span	10 mm
Span MIN	3 mm
Material clamp	Polypropylene (PP)
for clamping to	Beakers Pistons

