



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® Universal stand clamp Alu, clamping width 0-75 mm, 5 pcs./pack


Special Price

**€124.00** was

**€148.90**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Universal aluminum tripod clamp with cork insert

## Additional Information

---

No.	E-1395
Manufacturer (Brand)	neoLab
EAN	4058072064495
Transport temperature	Room temperature
Color	silver
Material	Aluminum
Form	pincer-shaped. round with 2x2 fingers
Length	225 mm
suitable for	Beakers Electrodes Centrifuge tubes
MAX span	75 mm
Span MIN	0 mm
for DM Tripod Rod MAX	13 mm
for DM tripod rod MIN	11 mm
for clamping to	Tripod rods

