



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Socorex stand for Calibrex™universal 525/530
for remote aspiration, dispenser capacity 100
ml**

**€247.00
plus VAT &
Shipping**

Product Images



Description

Holders for dispensers in various designs

Holder for remote and tube aspiration (Art. No. 7-7475 and 7-7476)

This dispenser holder facilitates aspiration from remote containers (distance < 10 m, height < 2 m).

Additional Information

No.	SR-0200
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7476
EAN	4058072367183
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
suitable for	Socorex Dispenser
Material plate	Steel
Material rod / arm	Steel

