



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex holder for Calibrex™ 525/530 dispensers, capacity 100 ml

€147.89
plus VAT &
Shipping

Product Images



Description

Holders for dispensers in various designs

Stand for increased stability (Art. No. 7-7474 and 7-7890)

Stability, especially practical for smaller bottles with insufficient stability.

Additional Information

No.	SR-0423
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7890
EAN	4058072369361
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
suitable for	Socorex Dispenser
Material plate	Steel
Material rod / arm	Steel

