



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® stand ring 120 mm ø

€9.60
plus VAT &
Shipping

Product Images



Description

The **ecoLab stand rings** made of powder-coated steel are robust and ideal for hanging separating funnels and other conical devices. Select the ring with the appropriate diameter for your needs. This allows you to separate phases without worry over an extended period of time.

- 12 mm stem diameter

Additional Information

No.	DT-0616
Manufacturer (Brand)	neoLab
old neoLab article no.	E-0135
EAN	4058072228927
Transport temperature	Room temperature
Color	silver
Material	Steel galvanized
Form	annular. closed
suitable for	Separating funnel
MAX span	120 mm
Span MIN	120 mm
for DM Tripod Rod MAX	12 mm
for DM tripod rod MIN	8 mm
for clamping to	Tripod rods

